INSTRUCTIONS FOR USE:
This package insert is not a reference to surgical techniques. It is designed to assist in using this product.

Product Description: **SURGISPON®**
(ABSORBABLE GELATIN SPONGE)

SURGISPON - Absorbable Haemostatic Gelatin Sponge is a sterile malleable, water-insoluble, gelatin absorbable sponge intended for haemostatic use by applying to a bleeding surface. SURGISPON is non-pyrogenic and biocompatible.

SURGISPON is a surgical haemostatic sponge, manufactured from highly purified first grade gelatin material for use in various surgical procedures. When implanted in vivo and used in appropriate amounts, it is completely absorbed within 3-4 weeks. When applied to bleeding mucosal regions, it liquefies within 2 to 5 days.

SURGISPON gelatin sponges have a porous structure which activates the thrombocytes at the moment blood comes in contact with the matrix of the sponge. This causes the thrombocytes to release a series of substances which promote their aggregation at the same time as their surfaces change character, thus enabling them to act as a catalyst for the formation of the fibrin.

INDICATIONS:

SURGISPON can be effectively used in various surgeries for haemostasis, when control of capillary, venous and arteriole bleeding by pressure, ligation and other conventional procedures is ineffective or impractical.

CONTRA-INDICATIONS:

SURGISPON should not be used:

- in closure of skin incisions because it may interfere with the healing of skin edges
- in patients with known allergies to collagen
- in intra vascular compartments because of the risk of embolization
- in conjunction with methyl methacrylate adhesives
- in conjunction with autologous blood salvage circuits
- for the primary treatment of coagulation disorders
- in the presence of infection
- for controlling post-partum bleeding or menorrhagia
- in instances of pumping arterial haemorrhage
- Where blood or other fluids have pooled or in cases where the point of haemorrhage is submerged

SURGISPON Sponge will not act as a tampon or plug in a bleeding site, nor will it close off an area of blood collecting behind a tampon. The safety and effectiveness of SURGISPON have not been established:

- for use in ophthalmic procedures
- in children and pregnant women

The safety and effectiveness for use in urological procedures have not been established through a randomized clinical study. In urological procedures, gelatin sponge should not be left in the renal pelvis, renal calyces, bladder, urethra or ureters to eliminate the potential foci for calculus formation. Although the safety and

effectiveness of the combined use of gelatin sponge with other agents such as topical thrombin, antibiotic solution or antibiotic powder have not been evaluated in controlled clinical trials, if in the physician's judgement, concurrent use of topical thrombin or other agents is medically advisable, the product literature for that agent should be consulted for complete prescribing information.

HANDLING INSTRUCTIONS / WARNINGS / PRECAUTIONS:

SURGISPON may be used dry or saturated with physiological salt solution. Cut to the desired size, a piece of SURGISPON, either dry or saturated with sterile, isotonic sodium chloride solution (sterile saline), can be applied with light pressure directly to the bleeding site. When applied dry, a single piece of SURGISPON should be manually applied to the bleeding site, and held in place with moderate pressure until hemostasis results.

When used with sterile saline, SURGISPON should be first immersed in the solution and then withdrawn, squeezed between gloved fingers to expel air bubbles, and then replaced in saline until needed. The SURGISPON sponge should promptly return to its original size, with slight expansion in thickness and shape in the solution. If it does not, it should be removed again and kneaded vigorously until all air is expelled and it does expand to its original size, with slight increases in thickness and shape when returned to the sterile saline. SURGISPON if used wet it may be blotted to dampness on gauze before application to the bleeding site. It should be held in place with moderate pressure, using a pledget of cotton or small gauze sponge until hemostasis results. Removal of the pledget or gauze is made easier by wetting it with a few drops of sterile saline, to prevent pulling up the SURGISPON which by then should enclose a firm clot. Use of suction applied over the pledget of cotton or gauze to draw blood into the SURGISPON is unnecessary, as SURGISPON will draw up sufficient blood by capillary action.

The first application of SURGISPON will usually control bleeding, but if not, additional applications may be made. For additional applications, fresh pieces should be used, prepared as described above. Use only the minimum amount of SURGISPON necessary to produce hemostasis. Once haemostasis is achieved any excess Gelatin Sponge should be carefully removed because of the possibility of dislodgement of the device or compression of other nearby anatomic structures.

When placed into cavities or closed tissue spaces, minimal preliminary compression is advised. The over packing of SURGISPON, should be avoided, since recovery to its initial volume may interfere with normal function and/or could cause possible or eventual compression necrosis of surrounding tissue and nerve damage. While packing a cavity for haemostasis sometimes surgically indicated, Gelatin Sponge should not be used in this manner unless excess product not needed to maintain haemostasis is removed. SURGISPON should be removed after usage and bleeding has stopped in radical cavities, laminectomy procedures, around or in proximity to foramina in bone, areas of bony confine, the spinal cord and /or the optic nerve and chiasma or closed tissue spaces with presence of bone. This might lead to unintended pressure on neighboring structures which may result in pain for the patient or might create the potential for nerve damage.

Since SURGISPON causes little more cellular reaction than does the blood clot, the wound may be closed over it. SURGISPON may be left in place when applied to mucosal surfaces until it liquefies.

SURGISPON is not intended as a substitute for meticulous surgical technique and the proper application of ligatures or other conventional procedures for haemostasis.

SURGISPON should not be used for the primary treatment of coagulation disorders.

SURGISPON should be used with caution in contaminated areas of the body. If signs of infection or abscess develop where SURGISPON Sponge has been positioned, reoperation may be necessary in order to remove the infected material and allow drainage. Users should be familiar with surgical procedures and techniques involving gelatin sponge before employing SURGISPON. Particular factors of each patient should be considered regarding healing process in vivo. Contaminated or infected wounds should be managed with acceptable surgical practice.

SURGISPON sterile sponge is packed in sterile blisters/envelops which guarantee sterility. Once the package is opened, contents are subject to contamination. It is recommended that SURGISPON should be used as soon as the package is opened and unused contents discarded. It should not be reused or resterilized as it may lose its physical characteristics as well as sterility. Open the blister/envelops pack by pulling its two loose ends away from each other—in such a way that the SURGISPON falls out undamaged, onto a sterile surface.

The product is intended for single use! Discard any unused SURGISPON remaining. Dispose of contaminated devices and packaging materials utilizing standard hospital procedures and universal precautions for bio-hazardous waste.

STERILISATION:

SURGISPON is sterilised by gamma radiation. Do not resterilise! Do not use if the package is opened or damaged. This device was designed, tested and manufactured for single patient use only. Reuse, use of the device with opened or damaged packaging, reprocessing and/or resterilisation of this device may lead to its failure and subsequent injury, illness or death of the patient and/or create the risk of contamination and patient infection, illness or death of the patient.

STORAGE:

The product should be stored in its original packaging and the outer package must therefore be closed immediately after use. The box should be stored in a clean, dry room at a temperature not more than 30°C.

Do not use after the expiry date!

PACKING:

Available in various shapes and sizes suitable for different kinds of surgeries.

SYMBOLS USED ON THE LABEL:

-  Warning
-  Temperature limit
-  Keep dry
-  Keep away from sunlight
-  Do not re-use/for single use only
-  Do not resterilize
-  Do not use if package is damaged
-  Manufacturer
-  Use by date - Expiry Date: Month & Year
-  See instruction for use
-  Manufacturing Date
-  Batch code
-  Catalogue Number
-  Sterilized using irradiation

CE
2460

EC REP
OBELIS S.A (www.obelis.net)
Bd. Général Wahis, 53
1030 Brussels, Belgium

Manufactured by:
AEGIS LIFESCIENCES PVT. LTD.
215/216, Mahagujrat Industrial Estate,
Phase -III PO: Changodar, Gam: Moraiya,
Tal: Sanand, Distt: Ahmedabad 38 2213, India.
Web: www.aegis-lifesciences.com
Email: info@aegis-lifesciences.com
Ph.:+91 9726903131, +91 9727103131

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